



Introduction

Background:

Activation of hypoxia signaling pathways has been linked to cancer cell-intrinsic immune evasion and resistance. IFN γ and TNF α are immune cell-derived effector cytokines critical for tumor immunity. However, the effect of these cytokines in the hypoxic tumor microenvironment is less well studied.

Aim of the study:

To investigate how hypoxia modulates the kinase signaling pathways activated by the cytokines IFN γ and TNF α in melanoma cells, and to identify potential mechanisms underlying immune evasion in the tumor microenvironment.

Method:

Murine B16F10 melanoma cells were treated with 40 ng/ml of IFN-gamma and TNF-alpha under Normoxia (21 % oxygen) and hypoxia (1% oxygen) conditions for 6 or 24 hours. Tyrosine and Serine/threonine kinase activity profiles were generated using PamGene's KinomePro technology.

KinomePro Coverage




Fig.2. Kinases covered by PamChip® microarrays

Kinases affected upon Hypoxia




Fig.3. Kinases activated by hypoxia treatment of B16F10 cells

Pathways enriched by IFN- γ & TNF- α treatment




Fig.5. (A) Top enriched pathways heatmap
(B) PI3K pathway protein-protein interaction enrichment

PamGene's KinomePro technology

With the advanced PamChip® technology, we can simultaneously measure the activity of over 380 kinases, providing a comprehensive view of kinase signalling pathways in a wide range of cell types and tissues.




Fig.1. Kinase activity profiling on PamChip® microarrays

Tyrosine and Serine/Threonine kinases affected upon IFN- γ and TNF- α treatment




Fig.4. Protein tyrosine kinases (*left*) and serine/threonine kinases (*right*) differentially up-(red) or down-(blue) regulated upon IFN γ and TNF- α treatment for 6 hours (*top*) or 24 hours (*bottom*) in Normoxia or Hypoxia treated cells. Size of circle depicts specificity of kinase change.

Kinome profiling enables deeper insights into complex hypoxia biology and opens new avenues for drug target discovery